

Overview all Advanced Instruments osmometers | for Biopharma

<p>All models support Pharmacopeia testing guidelines.</p>				
<p>WORKFLOW FEATURES</p>	<p>OsmoTECH® Single-Sample Micro-Osmometer</p>	<p>OsmoTECH® XT Single-Sample Micro-Osmometer</p>	<p>OsmoTECH® PRO Multi-Sample Micro-Osmometer</p>	<p>OsmoTECH® HT Automated Micro-Osmometer</p>
<p>Best for labs with daily sample processing needs of</p>	<p>1-10/day</p>	<p>1-10/day</p>	<p>>10/day</p>	<p>>20/day</p>
<p>Sample volume</p>	<p>20 µL</p>	<p>20 µL</p>	<p>30 µL</p>	<p>50 µL</p>
<p>Sample introduction</p>	<p>One-step direct</p>	<p>One-step direct</p>	<p>Sample cup</p>	<p>96 well plate</p>
<p>Operating range (mOsm/kg H₂O)</p>	<p>0 - 2,000</p>	<p>0 - 4,000</p>	<p>0 - 2,000</p>	<p>0 - 2,000</p>
<p>Test time (seconds)</p>	<p>90</p>	<p>Less than 150 in low range, less than 180 in high range</p>	<p>90</p>	<p>140</p>
<p>COMPLIANCE SUPPORT FEATURES</p>				
<p>Supports 21 CFR part 11 & EU Annex 11 compliance</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>Supports Pharmacopeia testing guidelines compliance</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>INSTRUMENT FEATURES</p>				
<p>Integrated barcode scanner</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>On-board printer</p>	<p>✓*</p>	<p>✓*</p>	<p>✓*</p>	<p>✓*</p>
<p>Multi-language touchscreen display</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>Username & password compatibility</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>Multi-sample capacity</p>			<p>✓</p>	<p>✓</p>
<p>Automated cleaning between samples</p>			<p>✓</p>	<p>✓</p>
<p>Automated sample pipetting</p>				<p>✓</p>
<p>Temperature control method</p>	<p>Dry cooling (No heat transfer fluid)</p>	<p>Dry cooling (No heat transfer fluid)</p>	<p>Dry cooling (No heat transfer fluid)</p>	<p>Dry cooling (No heat transfer fluid)</p>
<p>DATA MANAGEMENT</p>				
<p>TCP/IP</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>USB</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>Webserver for remote data viewing</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>
<p>OPC-UA, Network share</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓*</p>
<p>Database backup</p>	<p>✓</p>	<p>✓</p>	<p>✓</p>	<p>✓*</p>
<p>DTE RS-232, Barcode port</p>	<p>✓</p>			

*Sold separately