




Overview all Advanced Instruments osmometers | for Clinical

<p>All models support HIPAA, CAP and CLIA compliance.</p>			
<p>WORKFLOW FEATURES</p>	<p>Osmo1® Single-Sample Micro-Osmometer</p>	<p>OsmoPRO® Multi-Sample Micro-Osmometer</p>	<p>OsmoPRO® MAX Automated Osmometer</p>
<p>Best for labs with daily sample processing needs of:</p>	<p>1 - 10/day</p>	<p>>10/day</p>	<p>>1/day</p>
<p>Sample volume</p>	<p>20 µL</p>	<p>20 µL</p>	<p>150 µL</p>
<p>Sample introduction</p>	<p>One-step direct</p>	<p>Sample cup</p>	<p>One-step direct</p>
<p>Operating range (mOsm/kg H₂O)</p>	<p>0 - 2,000</p>	<p>0 - 2,000</p>	<p>0 - 2,000</p>
<p>Test time (seconds)</p>	<p>90</p>	<p>90</p>	<p>180</p>
<p>COMPLIANCE SUPPORT FEATURES</p>			
<p>Supports CAP & CLIA compliance</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>
<p>Supports HIPAA compliance</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>
<p>INSTRUMENT FEATURES</p>			
<p>Integrated barcode scanner</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>
<p>On-board printer</p>	<p>☑</p>	<p>☑</p>	<p>☑*</p>
<p>Multi-language touchscreen display</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>
<p>Username & password compatibility</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>
<p>Multi-sample capacity</p>		<p>☑</p>	<p>☑</p>
<p>Automated cleaning between samples</p>		<p>☑</p>	<p>☑</p>
<p>Automated sample pipetting</p>			<p>☑</p>
<p>Temperature control method</p>	<p>Dry cooling (No heat transfer fluid)</p>	<p>Dry cooling (No heat transfer fluid)</p>	<p>Dry cooling (No heat transfer fluid)</p>
<p>DATA MANAGEMENT FEATURES</p>			
<p>LIS Connectivity</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>
<p>USB</p>	<p>☑</p>	<p>☑</p>	<p>☑</p>

*Sold separately