

# Technical Data

		SCI-Tive	SCI-tive Dual (Symmetrical)	SCI-tive Dual (Asymmetrical L=left chamber, R=right chamber)
<b>External Dimensions [mm]</b>	Width	1660	2610	2860
	Depth	826	980	970
	Hight	1077	1131	1850-1900
<b>Internal Dimensions [mm]</b>	Width	1200	1170	1170 ( L ) / 1000 ( R )
	Depth	600	555	555 ( L ) / 800 ( R )
	Hight	620	620	5590 ( L ) / 690 ( R )
<b>Maximum Capacity</b>	T75 Flasks	592	varies	varies
<b>Working Capacity</b>	T75 Flasks	180	varies	varies
<b>Interlock Dimensions [mm]</b>	Width	270	270	270
	Depth	200	200	200
	Hight	210	210	210
<b>Interlock Capacity</b>	T75 Flasks	15	15	15
<b>Interlock Time Cycle</b>	Seconds	60	60	60
<b>Weight</b>	kg	230	355	355